## **Printing instructions**

## 1 Print the sizing guide



- → Print on U.S. Letter- or A4-sized paper.
  → Ensure you select "Actual size" and not "Fit" in the printer settings.
- $\rightarrow$  For best results, print in color.



## Measure your printout

ightarrow Place your ruler on the scale below and ensure they measure the same.

Л 

## Sizing guide instructions





Hold the guide so that it presses against the underside of your nose.





Choose the size that best fits your nose.

F&P and Evora are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip. REF 620938 REV A 2021-02 © 2021 Fisher & Paykel Healthcare Limited.

